

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.7 Revision Date 27.02.2025 Print Date 04.05.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : Hu Ghrelin Act (BULKPKG) ELISA COLDPKS

Product Number : EZGRA
Catalogue No. : 6D2892
Brand : Millipore

REACH No. :

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3

CHEMIKART

1.4 Emergency telephone

Emergency Phone # : 000 800 1007 141 (CHEMTREC)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Corrosive to Metals, (Category 1) H290: May be corrosive to metals.

Acute toxicity, (Category 4) H302: Harmful if swallowed.

Acute toxicity, (Category 4) H332: Harmful if inhaled.

Acute toxicity, (Category 3) H311: Toxic in contact with skin.

Serious eye damage, (Category H318: Causes serious eye damage.

1)

Skin sensitization, (Category 1) H317: May cause an allergic skin reaction.

Millipore- EZGRA Page 1 of 12

Specific target organ toxicity - repeated exposure, (Category 2),

Respiratory Tract

H373: May cause damage to organs through prolonged or repeated exposure.

Long-term (chronic) aquatic hazard, (Category 2)

H411: Toxic to aquatic life with long lasting

effects.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal Word Danger

Hazard Statements

H290 May be corrosive to metals.

H302 + H332 Harmful if swallowed or if inhaled.

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H373 May cause damage to organs (Respiratory Tract) through

prolonged or repeated exposure.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsina.

P314 Get medical advice/ attention if you feel unwell.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram

Signal Word Danger

Hazard Statements

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

Precautionary Statements

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/

doctor if you feel unwell.

Millipore- EZGRA Page 2 of 12

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

2.3 Other hazards - none

SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

No metal containers.

Tightly closed. Keep in a well-ventilated place. Keep locked up or in an area accessible only to qualified or authorized persons.

Storage stabilityRecommended storage temperature

2 - 8 °C

Storage class

Storage class (TRGS 510): 6.1C: Combustible, acute toxic Cat.3 / toxic compounds or compounds which causing chronic effects

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

Millipore- EZGRA Page 3 of 12

SECTION 8: Exposure controls/personal protection

Refer to component SDS

SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport information

14.1 UN number

ADR/RID: 3316 IMDG: 3316 IATA: 3316

14.2 UN proper shipping name

ADR/RID: CHEMICAL KIT IMDG: CHEMICAL KIT IATA: Chemical kit

14.3 Transport hazard class(es)

ADR/RID: 9 IMDG: 9 IATA: 9

14.4 Packaging group

ADR/RID: III IMDG: III IATA: III

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

Tunnel restriction code : (E)

Further information : No data available

Millipore- EZGRA Page 4 of 12

14.7 Maritime transport in bulk according to IMO instruments

Not applicable for product as supplied.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 16: Other information

Kit Components: Microtiter ELISA Plate 10X HRP Wash Buffer Human Ghrelin (Active) Standard lyophilized (2ng/ml) COLD PACKS	Millipore Millipore Millipore	EPRAM EWB-HRP E8088-K	Skin Sens. 1; H317 Acute Tox. 4; Acute Tox. 3; Eye Dam. 1; STOT RE 2; Aquatic Chronic 2; H302 + H332, H311, H318, H373, H411
Human Ghrelin (Active) Quality Controls 1,2 Iyophilized (0.5ml) COLD PACKS	Millipore	E6088-K	Acute Tox. 4; Acute Tox. 3; Eye Dam. 1; STOT RE 2; Aquatic Chronic 2; H302 + H332, H311, H318, H373, H411
Matrix Solution (1 ml) Assay Buffer (15 ml) COLD PACKS	Millipore Millipore	EMTX-GA EABGR	Met. Corr. 1; H290
Human Ghrelin (Active) Capture Antibody (3 ml) COLD PACKS	Millipore	E1088-C	Eye Irrit. 2; H319
Human Ghrelin (Active) Detection Antibody (3 ml) COLD PACKS	Millipore	E1088-D	Eye Irrit. 2; H319
Enzyme Solution Substrate (TMB) (12 ml)	Millipore Millipore	EHRP-88 ESS-TMB2	Skin Sens. 1; H317
Stop Solution (12 ml)	Millipore	ET-TMB-M	Met. Corr. 1; H290

Millipore- EZGRA Page 5 of 12

Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the	mixture	Classification procedure:
Met. Corr.1	H290	Calculation method
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.3	H311	Calculation method
Eye Dam.1	H318	Calculation method
Skin Sens.1	H317	Calculation method
STOT RE2	H373	Calculation method
Aquatic Chronic2	H411	Calculation method

Further information

The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact

Millipore- EZGRA Page 6 of 12

with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.



Millipore- EZGRA Page 7 of 12

Kit Components:

Microtiter EL	Microtiter ELISA Plate				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
EPRAM					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

10X HRP Was	h Buffer				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
identifier Millipore - EWB-HRP		Warning	May cause an allergic skin reaction.	Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water. Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical advice/ attention. Wear protective gloves. Take off contaminated clothing and wash it before	information
				reuse.	

Human Ghrelin (Active) Standard lyophilized (2ng/ml) COLD PACKS					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
E8088-K			swallowed or if	the	
			inhaled.	environment. IF	
			Toxic in contact	ON SKIN: Wash	
			with skin.	with plenty of	
			Causes serious	water.Call a	
			eye damage.	POISON	

Millipore- EZGRA Page 8 of 12

May cause	CENTER/ doctor
damage to	if you feel
organs	unwell. Wear
(Respiratory	protective
Tract) through	gloves/
prolonged or	protective
repeated	clothing/ eye
exposure.	protection/ face
Toxic to aquatic	protection, Idea
life with long	INHALED:
lasting effects.	Remove person
lasting chects.	to fresh air and
	keep
	comfortable for
	breathing. Call a
	POISON
	CENTER/ doctor
	if you feel
	unwell. IF IN
	EYES: Rinse
	cautiously with
	water for
	several minutes.
	Remove contact
	lenses, if
	· · · · · · · · · · · · · · · · · · ·
	present and easy to do.
	Continue
	rinsing. Get
	medical advice/
	attention if you
	feel unwell.

Human Ghreli	Human Ghrelin (Active) Quality Controls 1,2 lyophilized (0.5ml) COLD PACKS				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
E6088-K			swallowed or if	the	
			inhaled.	environment. IF	
			Toxic in contact	ON SKIN: Wash	
			with skin.	with plenty of	
			Causes serious	water.Call a	
			eye damage.	POISON	
			May cause	CENTER/ doctor	
			damage to	if you feel	
			organs	unwell. Wear	
			(Respiratory	protective	
			Tract) through	gloves/	
			prolonged or	protective	
			repeated	clothing/ eye	
			exposure.	protection/ face	
			Toxic to aquatic	protection. IF	
			life with long	INHALED:	
			lasting effects.	Remove person	
				to fresh air and	

Millipore- EZGRA Page 9 of 12

	keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/ attention if you	
	feel unwell.	

Matrix Solution (1 ml)					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Warning	May be	Keep only in	
EMTX-GA			corrosive to	original	
			metals.	packaging.	
				Absorb spillage	
				to prevent	
				material	
				damage.	

Assay Buffer	Assay Buffer (15 ml) COLD PACKS				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
EABGR					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

Human Ghreli	Human Ghrelin (Active) Capture Antibody (3 ml) COLD PACKS					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore -		Warning	Causes serious	Wash skin		
E1088-C			eye irritation.	thoroughly after		
				handling. IF IN		
				EYES: Rinse		
				cautiously with		
				water for		
				several minutes.		
				Remove contact		

Millipore- EZGRA Page 10 of 12

lenses, if	
present and	
easy to do.	
Continue	
rinsing. Wear	
eye protection/	
face protection.	
If eye irritation	
persists: Get	
medical advice/	
attention.	

Human Ghrelin (Active) Detection Antibody (3 ml) COLD PACKS					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Warning	Causes serious	Wash skin	
E1088-D			eye irritation.	thoroughly after	
				handling. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Wear	
				eye protection/	
				face protection.	
				If eye irritation	
				persists: Get	
				medical advice/	
				attention.	

Enzyme Solution					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Warning	May cause an	Avoid breathing	
EHRP-88			allergic skin	dust/ fume/	
			reaction.	gas/ mist/	
				vapors/ spray.	
				IF ON SKIN:	
				Wash with	
				plenty of water.	
				Contaminated	
				work clothing	
				should not be	
				allowed out of	
				the workplace.	
				If skin irritation	
				or rash occurs:	
				Get medical	
				advice/	
				attention.	

Millipore- EZGRA Page 11 of 12

	Wear protective gloves. Take off contaminated clothing and wash it before
	reuse.

Substrate (TMB) (12 ml)					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - ESS-TMB2					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Stop Solution (12 ml)					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - ET-		Warning	May be	Keep only in	
TMB-M			corrosive to	original	
			metals.	packaging.	
				Absorb spillage	
				to prevent	
				material	
				damage.	

Millipore- EZGRA Page 12 of 12